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Cell Viability & Cell Death Assay Tools

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
產品功能	Product	Detection Method
檢測粒線體脫氫酶 (Dehydrogenase) 活性	Cell Counting Kit 8 (WST-8 / CCK8) 熱銷推薦	Microplate Reader (Colorimetric)
	WST-1 Assay Reagent	Microplate Reader (Colorimetric)
	XTT Assay Kit	Microplate Reader (Colorimetric)
	MTS Assay Kit	Microplate Reader (Colorimetric)
	MTT Reagent	Microplate Reader (Colorimetric)
檢測細胞內 ATP 含量	ATPlite Assay Kit 熱銷推薦 (詳見下頁說明)	Microplate Reader (Luminescent)
	ATPlite 1step Assay Kit 熱銷推薦 (詳見下頁說明)	Microplate Reader (Luminescent)
	ATP Assay Kit	Microplate Reader (Colorimetric/Fluorometric)
	ADP/ATP Ratio Assay Kit	Microplate Reader (Luminescent)

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
產品功能	Product	Detection Method
檢測細胞膜破損 (Cell membrane damage)	LDH Assay Kit 	Microplate Reader (Colorimetric)
	DELFI [®] EuTDA Cytotoxicity Reagents 	Microplate Reader (TRF)
	Propidium Iodide	Fluorescence microscope/ Flow cytometry
	DRAQ7™	Fluorescence microscope/ Flow cytometry
	Trypan Blue Staining Solution	Light microscope

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
DELFI[®] EuTDA Cytotoxicity Reagents (continued)

Assay Principle of DELFIA[®] EuTDA Cytotoxicity Reagents. The DELFIA cytotoxicity assay takes advantage of a fluorescence enhancing ligand (BATDA) which crosses the cell membrane passively, enabling target cell loading that is both rapid and gentle. Once inside the cell, the ligand is immediately hydrolyzed by cellular esterases to generate a hydrophilic molecule (TDA) that can no longer penetrate the membrane. Cytolysis releases TDA into the supernatant to form a highly fluorescent lanthanide chelate with Europium; cytotoxicity levels are proportional to the amount of fluorescent signal produced.

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
產品功能	Product	Detection Method
檢測 DNA 複製 (DNA replication) 與細胞分裂 (Cell division) 活化程度	BrdU Cell Proliferation ELISA Kit 熱銷推薦	Microplate Reader (Colorimetric)
	EdU Assay / EdU Staining Proliferation Kit (iFluor 488)	Fluorescence microscope/ Flow cytometry
	Propidium Iodide Flow Cytometry Kit	Flow cytometry
	DRAQ5™	Fluorescence microscope/ Flow cytometry
	CFSE - Cell Labeling Kit	Fluorescence microscope/ Flow cytometry

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
Paraffin fixed human tonsil tissue sections (10 μ m) are stained with TUNEL Assay Kit - HRP-DAB. (A) Counter-stained with methyl green; (B) Eliminate counter-stain; (C) Section treated with DNase I as positive control; (D) Negative control without TdT enzyme.

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
產品功能	Product	Detection Method
檢測老化特異性蛋白 (Senescence-associated beta galactosidase)	Senescence Detection Kit  熱銷推薦	Light Microscope
	Senescence Assay Kit (Beta Galactosidase, Fluorescence)	Flow cytometry

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Senescence-associated beta-galactosidase staining at the (A) second, (B) fourth, (C) eighth, (D) twelfth passage *in vitro* expansion.

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