

# 免疫螢光染色 細胞影像觀察耗材 買2送1

限本頁產品 (可互搭)  
2019/2/18 - 3/26




## 固定樣本影像觀察



## Chamber Slides

- 由矽膠材質透明框架與玻片所組成，撕下框架輕鬆不費力。
- 可同時進行 3 / 8 / 12 組平行測試實驗。
- 3 well chamber 另提供專用玻片，可在節省抗體染劑用量 (150  $\mu$ l/well) 的同時仍保有大範圍觀察面積。
- 適合細胞培養與免疫螢光染色實驗使用，染色結果可以一般蓋玻片封片後長期保存。
- 適用於正立式及倒立式顯微鏡觀察。

## 活細胞影像觀察



## $\mu$ -Slide and $\mu$ -Dish

- 超薄底設計，厚度符合 No.1.5 標準蓋玻片厚度，上高倍率顯微鏡觀察也沒問題。
- ibiTreat 底部適合各類細胞貼附及生長。
- 極佳的物理及化學耐受性，相容於多種常見之固定試劑 (methanol, acetone, paraformaldehyde)。

### Chamber Slides

Product	Size	Cat. No.
3 Well Chamber, removable	15	IB-80381
8 Well Chamber, removable	15	IB-80841
12 well Chamber, removable	15	IB-81201

### $\mu$ -Slide

Product	Size	Cat. No.
$\mu$ -Slide 2 well ibiTreat	15	IB-80286
$\mu$ -Slide 2 well Ph+ ibiTreat	15	IB-80296
$\mu$ -Slide 4 well ibiTreat	15	IB-80426
$\mu$ -Slide 4 well Ph+ ibiTreat	15	IB-80446
$\mu$ -Slide 8 well ibiTreat	15	IB-80826
$\mu$ -Slide 2 x 9 well ibiTreat	15	IB-81806

### $\mu$ -Dish

Product	Size	Cat. No.
$\mu$ -Dish 35 mm, low ibiTreat	60	IB-80136
$\mu$ -Dish 35 mm, low ibiTreat Grid-500	60	IB-80156
$\mu$ -Dish 35 mm, high ibiTreat	60	IB-81156
$\mu$ -Dish 35 mm, high ibiTreat Grid-500	60	IB-81166
$\mu$ -Dish 35 mm, high Glass Bottom	60	IB-81158
$\mu$ -Dish 35 mm, high Glass Bottom Grid-50	30	IB-81148
$\mu$ -Dish 35 mm, high Glass Bottom Grid-500	30	IB-81168
$\mu$ -Dish 35 mm, high Bioinert	30	IB-81150
$\mu$ -Dish 35 mm, Quad ibiTreat	60	IB-80416

680+  
Publications  
about ibidi  
IF products

ibidi IF  
products  
overview

