

Data Sheet

Product Name: Multiple Species Paraffin Tissue Array - Liver

Catalog No.: T8134149

Lot No.: A805189

Species: Human Mouse Rat Monkey Plant Cell line Multiple Species

Tissue Type: Normal Adult Fetal Tumor Disease

Tissue Array Diagram:

	1	2	3	4	
A	x	O	O	O	
B	O	O	O	O	
C	O	O	O	O	
D	O	O	O	O	

Donor Information

Position	Species	Race/Strain/Breed	Tissue Type	Age	Sex	Diagnosis
A1	Blank					
A2	Human	Caucasian	Liver	70 years	F	Normal
A3	Mouse	Balb/C	Liver	8.5 weeks	F	Normal
A4	Rat	Wistar	Liver	4.5 weeks	F	Normal
B1	Monkey	Rhesus	Liver	1 year	F	Normal
B2	Monkey	Cynomolgus	Liver	3.2 year	F	Normal
B3	Rabbit	New Zealand White	Liver	7.5 weeks	F	Normal
B4	Guinea pig	guinea pig	Liver	5.5 weeks	F	Normal
C1	Goat	Bora	Liver	7month	F	Normal
C2	Sheep	Dorset	Liver	7month	F	Normal
C3	Hamster	Golden Syrian Hamster	Liver	8 weeks	F	Normal
C4	Chicken	Leghorn Hens	Liver	Adult	F	Normal
D1	Bovine (cow)	Cross breed	Liver	2 year	F	Normal
D2	Horse	Cross breed	Liver	Adult	F	Normal
D3	Dog	Beagle	Liver	1.5 year	F	Normal
D4	Porcine (pig)	Cross breed	Liver	7.5 month	F	Normal

Components:

- 2 Tissue Array Slides 5 Tissue Array Slides per package
- Certificate of Analysis
- Data Sheet

FOR RESEARCH USE ONLY

APPROVED BY: _____