

Data Sheet

Product Name: Tissue Array

Catalog No.: T8234708-5

Lot No.: B101154

Species: Human Mouse Rat Monkey Plant Cell line

Tissue Type: Normal Adult Fetal Tumor Disease

Tissue Array Diagram:

	A	B	C	D	E	F	G	H	I	J	K	
1	O	O	O	O	O	O	O	O	O	O	O	
2	O	O	O	O	O	O	O	O	O	O	O	
3	O	O	O	O	O	O	O	O	O	O	O	
4	O	O	O	O	O	O	O	O	O	O	O	
5	O	O	O	O	O	O	O	O	O	O		
6	O	O	O	O	O	O	O	O	O	O		

Components:

1. 5 Tissue Array Slides per package
2. Certificate of Analysis

FOR IN VITRO RESEARCH USE ONLY

APPROVED BY: _____



Data Sheet

Product Name: Tissue Array
Catalog No.: T8234708-5
Lot No.: B101154
Donor Information

Position	Tissue Type	Age	Gender	Diagnosis	Tissue ID
A1, A2	Placenta	29	F	Normal	HT01155
B1, B2	Blank				
C1, C2	Adipose	26	M	Normal	HT01156
D1, D2	Bladder	36	F	Normal	HT01157
E1, E2	Brain	68	M	Normal	HT01158
F1, F2	Breast	27	F	Normal	HT01159
G1, G2	Cerebellum	73	M	Normal	HT01160
H1, H2	Cervix	30	F	Normal	HT01161
I1, I2	Colon	28	M	Normal	HT01162
J1, J2	Diaphragm	26	M	Normal	HT01163
K1, K2	Duodenum	24	M	Normal	HT01164
A3, A4	Esophagus	23	M	Normal	HT01165
B3, B4	Gallbladder	24	M	Normal	HT01166
C3, C4	Heart	27	M	Normal	HT01167
D3, D4	Ileum	62	M	Normal	HT01168
E3, E4	Jejunum	26	M	Normal	HT01169
F3, F4	Kidney	56	M	Normal	HT01170
G3, G4	Liver	30	M	Normal	HT01171
H3, H4	Lung	24	M	Normal	HT01172
I3, I4	Ovary	50	F	Normal	HT01173
J3, J4	Pancreas	50	M	Normal	HT01174
K3, K4	Thyroid	26	M	Normal	HT01175
A5, A6	Rectum	24	F	Normal	HT01176
B5, B6	Skin	83	F	Normal	HT01177
C5, C6	Skeletal Muscle	26	M	Normal	HT01178
D5, D6	Spleen	30	M	Normal	HT01179
E5, E6	Stomach	66	M	Normal	HT01180
F5, F6	Testis	32	M	Normal	HT01181
G5, G6	Thymus	28	M	Normal	HT01182
H5, H6	Smooth Muscle	23	M	Normal	HT01183
I5, I6	Tonsil	28	F	Normal	HT01184
J5, J6	Uterus	40	F	Normal	HT01185