

# 秋季抗體特惠

Abcam 一級抗體買就送免費贈品 (WB 正貨產品擇一), 或享壹千元超值加購優惠 (WB 明星商品擇一)! 活動時間: 2020年10月15日 - 11月25日



## 免費贈品區 (任選一)

Type	Cat. No.	Product
Protease Inhibitor Cocktail	ab270055	Protease Inhibitor Cocktail EDTA Free, 3ml
	ab271306	Protease Inhibitor Cocktail, 3ml
Buffer	ab126587	10X Blocking Buffer, 50ml
	ab156034	10X RIPA Buffer, 15ml
	ab270550	Western Blot Stripping Buffer, 500ml (詳見背頁說明)
Protein Ladder	ab116028	Prism Ultra Protein Ladder (10-245 kDa), 500µl
	ab116029	Prism Ultra Protein Ladder (3.5 - 245 kDa), 500µl
	ab116027	Prism Ultra Protein Ladder (10 - 180 kDa), 500µl
Transfer Membrane	ab133411	Low Fluorescence PVDF Membrane (0.45µm, 7x9cm), 10 sheets
	ab133412	Nitrocellulose Transfer Membrane (0.45µm, 7x9cm), 10 sheets
	ab133413	Nitrocellulose Transfer Membrane (0.22µm, 7x9cm), 10 sheets
Staining Solution	ab270042	Ponceau S Solution, 1000ml
ECL Substrate	ab133408	ECL Substrate Kit (Very High Sensitivity), 20ml
	ab133409	ECL Substrate Kit (Ultra High Sensitivity), 20ml
Cell Lysates, Fractions, Extracts	ab136806	HeLa Apoptosis Lysate Set: Staurosporine-Treated and Vehicle-Treated Control, 2units
	ab136811	HeLa Cell Cycle Lysates: Thymidine-Treated, Nocodazole-Treated & Untreated Asynchronous Control, 3units
	ab168542	HeLa Cell Fractions (Cytoplasmic/Membrane/Nuclear), 3units
	ab157396	HeLa DNA Damage Whole Cell Lysate Set: UV Treated and Untreated Control, 2units
	ab139414	HeLa Histone Acetylation-Inhibited Whole Cell Lysate Set: Trichostatin A-Treated and Vehicle-Treated, 2units
	ab154479	HeLa Mitochondrial DNA Depletion Cell Lysate Set: Ethidium Bromide-Treated and Vehicle-Treated, 2units
	ab135395	HeLa Ser/Thr Phosphatase-Inhibited Whole Cell Lysate Set: Calyculin-Treated and Vehicle-Treated, 2units
	ab170195	HeLa Signaling Cascade Whole Cell Lysate Set: PMA-Treated and Vehicle-Treated Control, 2units
	ab170197	HeLa Whole Cell RIPA Extract, 2units



## 壹千元超值加購區 (任選一)

Type	Cat. No.	Product
Precast Gel*	ab270467	Bis-Tris Precast Gels - RunBlue™ (4-12%, 12-well, 10x10cm), 10units*
	ab270469	Bis-Tris Precast Gels - RunBlue™ (4-12%, 17-well, 10x10cm), 10units*
ECL Substrate	ab270531	ECL Pico - LumiBlue™, 250ml
	ab270542	ECL Express - LumiBlue™, 250ml
Bradford Reagent	ab119216	Bradford Reagent, 500ml (詳見背頁說明)

\*加購【Precast Gel】明星產品, 即可免費獲得加碼贈品【Running Buffer】(以下兩款 Buffer 產品任選一)

- MES Run Buffer (20X) - RunBlue™ (Bis-Tris), 500ml (ab270224)
- MOPS Run Buffer (20X) - RunBlue™ (Bis-Tris), 500ml (ab270225)



連續五年獲得 CiteAb 獎項肯定  
<https://www.citeab.com/awards>



Knockout (KO) 驗證抗體  
提供更準確可靠的實驗結果



RabMAb® 專利兔單株抗體技術  
提供更優越穩定的實驗再現性




63,000+ Primary Antibodies



# Western Blot Stripping Buffer

Cat. No.: ab270550 | Pack Size: 500 ml, 5 L

- 不需加熱步驟，室溫 20 分鐘即可迅速完成抗體脫附！
- 效果穩定，不會造成轉漬膜上的抗原損傷。
- 可進行多次抗體探測，有效節省珍貴樣本。



**Saves sample by re-probing.** Reprobing of the same membrane with different antibodies after treatment with Western Blot Stripping Buffer (ab270550).




# Bradford Reagent

Cat. No.: ab119216 | Pack Size: 500 ml

- 專利改良配方，定量結果更準確，不易受到樣本中的 Detergent 影響！  
※高含量定量流程，可耐受 1% Detergent；微量定量流程，可耐受 0.1% Detergent。
- 偵測範圍：0.1~1.5 mg/ml (高含量定量) 或 1~25 µg/ml (微量定量)。

Leading Competitor Bradford Assay

Abcam Bradford Reagent Assay



Abcam Bradford Reagent is not affected by the detergents compared to leading competitor "Coomassie Protein Assay" which are subject to the high background noise at 595 nm in the presence of detergents.

Western Blot Stripping Buffer 與 Bradford Reagent 皆為 Expecon (part of Abcam) 所開發的高效蛋白質實驗分析試劑。完整產品資訊與最新活動訊息，請洽 Abcam 台灣代理 — 伯森生技。



## 相關文獻：

**Bradford Reagent (ab119216)**  
previously called Bradford Ultra (BFU05L, Expecon)

- PLRP2 selectively localizes synaptic membrane proteins via acyl-chain remodeling of phospholipids. *J Lipid Res.* 2020 Sep 22;61(9):2833-2843. PMID: 32963038
- Bacterial polyphosphates interfere with the innate host defense to infection. *Nat Commun.* 2020 Aug 12;11(1):4035. PMID: 32788578
- The RNA-Binding Protein A1CF Regulates Hepatic Fructose and Glycerol Metabolism via Alternative RNA Splicing. *Cell Rep.* 2019 Oct 8;29(2):283-300.e8. PMID: 31597092
- Aquaporins influence seed dormancy and germination in response to stress. *Plant Cell Environ.* 2019 Aug;42(8):2325-2339. PMID: 30986891

## 相關連結：



Western Blot Stripping Buffer (ab270550)



Bradford Reagent (ab119216)

