


實驗結果難以重複？

如同《Nature》指出的，目前生科研究領域正面臨著實驗結果難以重複再現的嚴峻挑戰，而【**重組抗體**】即被視為解決此難題的最佳解答 [PMID [25652980](#)]。

Abcam 目前擁有近 19,000 項**重組抗體產品**，以精確快速的抗體製程技術與高規嚴謹的品質確校標準，為您帶來更多核心優勢，協助您更快更好地達成研究目標！

- **更高的一致性和再現性** —— 重組抗體是由一組特定的基因序列發展而來，所以生產過程是可控且可靠的，能夠避免一般單株抗體融合瘤細胞製程中的問題，如基因流失、基因突變和細胞株漂變等。這使得重組抗體批次間的差異非常小，因此能夠提供高度可重複的實驗結果。
- **更高的靈敏度和專一性** —— 透過重組技術，更容易使用基因工程提高抗體的專一性和靈敏度，並且能在選殖過程中篩選出最佳的抗體。
- **配方多元** —— 特別推出不含 BSA、Azide、Glycerol 等物質的無載體 (carrier-free) 配方抗體，使抗體能達到最佳的標定效率。
- **多種宿主動物選擇** —— 除兔子以外，尚有小鼠、大鼠、人類、山羊、羊駝、倉鼠等。




使用不同批號的 PD-L1 重組抗體進行 IHC 偵測，實驗結果展現高度再現性。

Recombinant Anti-PD-L1 antibody [28-8] (ab205921) References (137)



重組抗體比一般單株抗體更加高效穩定，並具有多種宿主動物選擇。



抗體買 1 送 1 即日起訂購一級抗體送重組抗體贈品

❖ 一級抗體金額需滿 15600 元方可參加本活動 ❖ 重組抗體贈品請見本刊背頁 ❖ 活動至 2020 年 2 月 27 日截止 ❖ 活動詳情請洽 Abcam 台灣代理一伯森生技



重組抗體贈品區

Cat. No.	Recombinant Antibody Product Name	Host	Size	Application
ab213204	Anti- 6X His tag ® antibody [EPR20547] - ChIP Grade	Rabbit	10 µl	ChIP, Flow Cyt, ICC/IF, IHC-P, IP, WB
ab176560	Anti- alpha Tubulin antibody [EPR13478(B)] - Loading Control	Rabbit	10 µl	Flow Cyt, ICC/IF, IHC-P, IP, WB
ab124979	Anti- Cofilin antibody [EPR6375] - Loading Control	Rabbit	10 µl	ICC/IF, IHC-P, WB
ab202554	Anti- COX IV antibody [EPR9442(ABC)] - Mitochondrial Loading Control	Rabbit	10 µl	Flow Cyt, ICC/IF, IHC-P, IP, WB
ab178397	Anti- Cyclophilin B antibody [EPR12703(B)] - Loading Control	Rabbit	10 µl	WB
ab205606	Anti- DDDDK tag (Binds to FLAG® tag sequence) antibody [EPR20018-251]	Rabbit	10 µl	Flow Cyt, ICC/IF, IHC-P, IP, WB
ab181602	Anti- GAPDH antibody [EPR16891] - Loading Control	Rabbit	10 µl	Flow Cyt, ICC/IF, IHC-P, IP, WB
ab183734	Anti- GFP antibody [EPR14104]	Rabbit	10 µl	Flow Cyt, ICC/IF, IHC-P, WB
ab111947	Anti- GST antibody [EPR4236]	Rabbit	10 µl	ICC/IF, IP, WB
ab182009	Anti- HA tag antibody [EPR4095]	Rabbit	10 µl	WB
ab45134	Anti- Hsp60 antibody [EP1006Y] - Loading Control	Rabbit	10 µl	IHC-P, WB
ab213511	Anti- mCherry antibody [EPR20579]	Rabbit	10 µl	ELISA, Flow Cyt, IP, WB
ab76020	Anti- Sodium Potassium ATPase antibody [EP1845Y] - Plasma Membrane Loading Control	Rabbit	10 µl	Flow Cyt, ICC/IF, IHC-P, WB
ab109503	Anti- Transferrin antibody [EPR2932(2)] - Serum Loading Control	Rabbit	10 µl	IHC-P, WB
ab154856	Anti- VDAC1 / Porin antibody [EPR10852(B)] - Mitochondrial Loading Control	Rabbit	10 µl	ICC/IF, IHC-Fr, IHC-P, WB
ab183497	Anti- VSV-G tag antibody [EPR12997]	Rabbit	10 µl	WB
ab109237	Anti- YY1 antibody [EPR4652] - Nuclear Loading Control and ChIP Grade	Rabbit	10 µl	ChIP, Flow Cyt, ICC/IF, IHC-P, IP, WB



抗體買1送1

即日起訂購一級抗體送重組抗體贈品

❖ 一級抗體金額需滿 15600 元方可參加本活動 ❖ 重組抗體贈品請見本頁說明 ❖ 活動至 2020 年 2 月 27 日截止 ❖ 活動詳情請洽 Abcam 台灣代理—伯森生技

